



# User Information Sheet

Ronco Item No : **229, 239, 249, 259**

Brand Name :

## RONCO PURE-TOUCH (LP)

Description : Synthetic Stretch Examination Gloves  
 Material : Blended of premium synthetic materials  
 Donning: Lightly Powdered  
 Characteristics: Disposable, Non-Sterile, Ambidextrous  
 Colour: White  
 Texture: Smooth



Medical  
Device  
License



Specification:

Size :	S	M	L	XL
Item Number :	<b>229</b>	<b>239</b>	<b>249</b>	<b>259</b>
Length : (mm) ± 10 mm	240	240	240	240
Width : (mm) ± 5 mm	85	95	105	115
Thickness :	3 mil			
Physical Properties :				
Tensile strength : minimum	10 MPa			
Elongation : minimum	400%			
Quality Level :				
Pinholes :	AQL 1.5			
Major Defects :	AQL 2.5			
Minor Defects :	AQL 4.0			

Packaging :

C1M = 100 gloves x 10 dispensers = 1,000 gloves per case

Product Features :

- Premium quality synthetic materials
- Superior elasticity
- Ease of finger movement
- Excellent strength
- Reduces hand fatigue
- Latex Free
- Non-sterile
- Beaded cuffs

Applications :

- Food processing
- Food service
- Hospitality
- Meat packing
- Dental
- Medical
- Pharmaceutical
- Veterinary
- Automotive
- Beauty / Aesthetics
- Electronics
- Environmental
- Industrial
- Laboratory
- Printing

Storage & Handling :

Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F ( 37°C )

Special Notes :

Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).