



User Information Sheet

Ronco Item No : **771, 773, 775, 777, 779, 781, 785**

Brand Name : **RONCO**

Product Name : **COVA-CAP**

Description : Pleated Bouffant Caps

Material : Non-woven Polypropylene

Size: 21" / 24"

Colour: White / Yellow / Blue / Pink / Green

Characteristics: Pleated
Latex-Free elastic double stitched



Specification :

Item Number :	771	773	775	777	779	781	785
Colour :	White	Yellow	Blue	Pink	Green	White	Blue
Length fully expanded :	530 ± 5 mm (21")					600 ± 5 mm (24")	
Length at rest :	210 ± 10 mm (8½")					240 ± 10 mm (9½")	
Width of pleat :	22 mm						

Packaging : C1M = 250 caps x 4 dispensers = 1,000 caps per master pack

- Product Features :
- 100% non-woven polypropylene
 - Pleated
 - Latex-Free
 - Soft elastic headband for maximum protection and hair retention
 - Fluid resistant
 - Lightweight and breathable
 - Ideal for department colour coding

- Applications :
- Food Processing
 - Food Service
 - Meat Packing
 - Fisheries/Seafood
 - Bakery
 - Beauty/Aesthetics
 - Medical facilities
 - Laboratory
 - Assembly line work

Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture and excessive heat above 100°F (37°C).

Special Notes :

Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).