



User Information Sheet

Ronco Item No : **6115, 6155**

Brand Name : **MASkin**

Description : N95 PARTICULATE RESPIRATOR

Material : 3-Ply PP spunbond outer layer, PP meltblown high filtration layer, PP needle punched nonwoven inner layer

Colour : White

Characteristics : 95% filtration effective against particulate aerosols free of oils



Specification :

Item number :	6115	6155
NIOSH approval number	TC#84A-5530	TC#84A-6766
Feature:	Regular	Valved
Length : (mm) ± 5 mm	110 mm	
Width : (mm) ± 5 mm	125 mm	
Elastic Head Band x 2	Latex-Free elastic - width 7 mm	
Nose Clip	Material	Aluminum
	Length	90 ± 10 mm

Packaging #6115 : 20 masks x 12 boxes = 240 masks / case
 #6155 : 10 masks x 12 boxes = 120 masks / case

Product Features :

- 95% filtration effective against particulate aerosols free of oils
- Conforms to NIOSH 42 CFR standard
- Low breathing resistance
- Aluminum nose-clip
- Soft Foam Nosepiece
- Latex-Free elastic ultrasonic welded to cup
- Adjusts to all facial contours
- Valved or Regular

Application :

- Medical / Examination
- Industrial / Manufacturing
- Carpentry / Drywall
- Food Handling
- Demolition / Renovations
- General Cleaning / Maintenance
- Laboratory
- Janitorial / Sanitation
- Waste Management

Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture and excessive heat above 100°F (37°C).

Special Notes :

Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).

