



ISO 9001:2000

“Tough gloves for tough jobs”



Nitrile disposable gloves

Beaded cuff

8 mil thick

12” long



Medical Device License

RONCO BLURITE

Latex disposable gloves

Beaded cuff

13 mil thick

12” long



Medical Device License

RONCO SILKTEX

EXTRA PROTECTION LONG

For more information or to receive a free sample, call RONCO now\*

Call RONCO today for information on your local authorized distributor



Your Guarantee of CLEAN & SAFE™ Products

Tel: (905) 660-6700
Fax: (905) 660-6903
Email: ronco@ronco.ca
www.ronco.ca



ISO 9001:2000