



# User Information Sheet

Ronco Item No : **819, 829, 839, 849, 859**

Brand Name :

**RONCO SILKTEX (PF)**

Description : Latex Examination Glove  
 Material : Natural Rubber Latex  
 Donning: Powder Free  
 Characteristics: Disposable, Non-Sterile, Ambidextrous  
 Colour: Natural  
 Texture: Fully Textured



Specification:

Size :		XS	S	M	L	XL
Item Number :		<b>819</b>	<b>829</b>	<b>839</b>	<b>849</b>	<b>859</b>
Length : (mm)	± 10 mm	240	240	240	240	240
Width : (mm)	± 5 mm	75	85	95	105	115
Thickness :		5 mil				
Physical Properties :						
	Tensile strength : minimum	21 MPa				
	Elongation : minimum	700%				
Pinhole AQL:		AQL 1.5				

Packaging : C1M = 100 gloves x 10 dispensers = 1,000 gloves per case

- Product Features :
- Premium quality latex
  - Excellent flexibility
  - First rate elasticity
  - Ease of finger movement
  - Reduces hand fatigue

Applications : Made from premium quality latex, Ronco Silktex is an excellent medical general purpose glove. Its high elasticity and ease of finger movement make its comfortable to wear. Its strength and latex formulation gives it its top performance. Applications include;

- Food processing
- Food service
- Hospitality
- Meat packing
- Dental care
- Medical examination
- Pharmaceutical
- Veterinary
- Beauty / Aesthetics
- Electronics
- Industrial
- Laboratory

Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F ( 37 °C ).

Caution: This product contains natural rubber latex that may cause allergic reactions in some users.

Special Notes : Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).